Degudent® U

Dental PFM alloy according to DIN EN ISO 22674 and 9693 for conventional ceramic materials for the manufacture of dental prostheses

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CE 0483

Colour: white Type: 4

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
Indications: Inlays, crowns, bridges of all spans, resin or normal expanding ceramic veneers

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
Alloy Composition %

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Warning! May damage the organs if exposed for long periods or repeatedly. Do not breathe dust/fume/gas/mist/vapour/aerosol. If you feel unwell, seek medical advice/attention. Dangerous component: indium

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UFI: HAH0-A0WH-800F-HY9R

Technical specifications Degudent® U

Density (g/cm ³)	18,0
Melting range °C	1260-1150
Hardness (HV 5/30) s,h,c/f	200, 240, 190/200
Pre-heat temperature °C	850
Casting temperature °C	1410
Additional heating time after melting (s) electrical 180, inductive 5-10, flame 5-10	
0,2% Elastic limit (MPa) s, h, c	470, 580, 490/490
Ultimate elongation (%) s, h, c	7, 3, 11/11
E-modulus (GPa)	105
Flux	p. ex. Pasta Flux
CTE 25-500 (600) °C	13,8 (14,0) µm/mK
Highest recommended firing temperature °C	1050
Soldering recommendation:	Orplid® Keramik Lot 1050 Orplid® Lot 760
Welding rod:	Degudent® U, 0,35 mm
Soft annealing: 950 °C / 15 min, quench in water	
Hardening: 600 °C / 15 min	

s = soft, h = hard, c/f = after casting/firing

Operating Instruction

Applying of casting sprues:

according to dental standards (see general operating instruction)

Investing:

phosphate bonded investments

Pre-heating:

850 °C, 30 to 120 min at end temperature, corresponding to muffle size

Casting:

1370-1420 °C, If necessary, increase by 20-30 °C with resistance-heated casting units, graphite crucible or with open flame with ceramic crucible.

Trimming:

With carbide milling and/or ceramical bonded abrasives. Always sand in one direction! Blasting at a flat angle with 110 µm aluminium oxide, pressure max. 2 bar.

Cleaning:

ultrasonic cleaning or steaming with steam device

Oxide firing:

5 min at 930 °C without vacuum

Safety information

Metal dust is hazardous to health. Use an air exhaustor system when trimming or separating. In rare cases, persons may be allergic on the ingredients of the alloy or experience discomfort as a result of an electrochemical reaction.

Interactions

In rare cases, the occlusal or interdental contact between different alloys may cause discomfort due to an electrochemical reaction.

Warranty

Our application-specific recommendations, irrespective of whether they are given vocal, in writing or practical instruction, are based on our own experience and tests and can only be recommendations. The products are subject to continuous improvement. We therefore reserve the right to make changes in design or composition. Please keep dry and store at room temperature.

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C. HAFNER GmbH + Co. KG