## Orplid® Keramik 4

Dental PFM alloy according to DIN EN ISO 22674 and 9693 for conventional ceramic materials for the manufacture of dental prostheses

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Colour: white **Type: 4**

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
Barva : bela **Type: 4**

**Indications:** Inlays, crowns, bridges of all spans, resin or normal expanding ceramic veneers

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
**Alloy Composition %**

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**Sestava zlitine %**

Warning! May damage the organs if exposed for long periods or repeatedly. Do not breathe dust/fume/gas/mist/vapour/aerosol. If you feel unwell, seek medical advice/attention. Dangerous component: indium

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s = soft, h = hard, c/f = after casting/firing

## Operating Instruction

## Applying of casting sprues:

according to dental standards (see general operating instruction)

## Investing:

phosphate bonded investments

## Pre-heating:

850 °C, 30 to 120 min at end temperature, corresponding to muffle size

## Casting:

1440 °C, If necessary, increase by 20-30 °C with resistance-heated casting units, graphite crucible or with open flame with ceramic crucible.

## Trimming:

With carbide milling and/or ceramical bonded abrasives. Always sand in one direction! Avoid high speeds and strong contact pressure (flaking). Blasting at a flat angle with 50 -125 µm aluminium oxide, pressure max. 2 bar.

## Cleaning:

ultrasonic cleaning or steaming with steam device

## Oxide firing:

5 min at 960 °C without vacuum

## Safety information

Metal dust is hazardous to health. Use an air exhaustor system when trimming or separating. In rare cases, persons may be allergic on the ingredients of the alloy or experience discomfort as a result of an electrochemical reaction.

## Interactions

In rare cases, the occlusal or interdental contact between different alloys may cause discomfort due to an electrochemical reaction.

## Warranty

Our application-specific recommendations, irrespective of whether they are given vocal, in writing or practical instruction, are based on our own experience and tests and can only be recommendations. The products are subject to continuous improvement. We therefore reserve the right to make changes in design or composition. Please keep dry and store at room temperature.

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C. HAFNER GmbH + Co. KG