Orplid® Keramik PF

Dental PFM alloy according to DIN EN ISO 22674 and 9693 for conventional ceramic materials for the manufacture of dental prostheses

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
Barva: svetlo rumena

Type: 4

Indications: Inlays, crowns, bridges of all spans, resin or normal expanding ceramic veneers

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
s = soft, h = hard, c/f = after casting/firing

Operating Instruction

Applying of casting sprues:

according to dental standards (see general operating instruction)

Investing:

phosphate bonded investments

Pre-heating:

800 °C, 30 to 120 min at end temperature, corresponding to muffle size

Casting:

1320 °C, If necessary, increase by 20-30 °C with resistance-heated casting units, graphite crucible or with open flame with ceramic crucible.

Trimming:

With carbide milling and/or ceramical bonded abrasives. Always sand in one direction! Blasting at a flat angle with 50 - 125 µm aluminium oxide, pressure max. 2 bar.

Cleaning:

ultrasonic cleaning or steaming with steam device

Oxide firing:

5 min at 940 °C without vacuum. It is imperative to remove the oxides. Preferably, acidify for 10 - 15 minutes in a warm and clean neacid bath or carefully blast with 50 - 110 µm aluminium oxide. Finally, steam off thoroughly.

Safety information

Metal dust is hazardous to health. Use an air exhaustor system when trimming or separating. In rare cases, persons may be allergic on the ingredients of the alloy or experience discomfort as a result of an electrochemical reaction.

Interactions

In rare cases, the occlusal or interdental contact between different alloys may cause discomfort due to an electrochemical reaction.

Warranty

Our application-specific recommendations, irrespective of whether they are given vocal, in writing or practical instruction, are based on our own experience and tests and can only be recommendations. The products are subject to continuous improvement. We therefore reserve the right to make changes in design or composition. Please keep dry and store at room temperature.

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C. HAFNER GmbH + Co. KG