EN Instructions for Use

- Self-curing compomer cement for the permanent cementation of AGC secondary parts in tertiary structures

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English

Description

AGC® Cem is a self-curing compomer cement for the permanent cementation of AGC secondary parts (telescope crowns) and AGC bar sleeves in tertiary structures (model casting). AGC Cem is suitable for intraoral and extraoral use.

Shade

yellow

Processing time

The processing and curing times depend on the ambient temperature. Once AGC Cem has been dispensed from the automix syringe, the following times apply:

	Room temperature approx. 23 °C	Intraoral
Processing time	Approx. 90 s	Approx. 60 s
Curing time (including processing time)	approx. 15 min	approx. 6 min

Note

The restoration must not be put under stress before the hardening process is complete.

Mixing ratio

AGC Cem is extruded from the automix syringe in an optimum 1:1 mixing ratio.

Composition

Ionomer glass in a matrix composed of various dimethacrylates, activator, catalyst and additives.
Filler content: 68 wt% % = 47 Vol. % (0.02–20 µm).

Indication

Intraoral and extraoral cementation in tertiary structures of

- AGC telescope crowns
- AGC bar sleeves

Contraindications

- Cementation of restoration on titanium bonding bases
- Use in combination with materials containing eugenol
- Any indications not stipulated in the working procedures

Side effects

No systematic side effects are known to date. In individual cases, allergic reactions to individual components have been reported. AGC Cem should not be used, however, in the case of known allergies to any of the constituents.

Interactions

Phenolic substances (e.g. wintergreen oil, eugenol) inhibit polymerization. Consequently, the application of materials containing these substances is to be avoided.

Application

- Carefully blast the surface of the AGC structure with 50 μ Al₂O₃ at 1 bar pressure.
- Also blast the tertiary structure (observe the manufacturer's instructions).
- The structure must be free from grease and dust after blasting. After cleaning, the bonding surface must not be contaminated under any circumstances.
- Apply a suitable bonding agent (e.g. Monobond® Plus) on the blasted structure surface (observe the manufacturer's instructions).

Note:

Do not use bonding systems that form a thick layer on the structure.

- AGC Cem is directly applied from the automix syringe using a mixing tip.
- After opening and removing the cap, push an automix mixing tip onto the syringe. Line up the notch on the tip with the corresponding notch on the syringe and turn it 90° until it locks into position.
- Apply a sufficient amount of AGC Cem directly from the mixing tip on all bonding surfaces of the tertiary structure.
- Press the parts lightly and evenly together and check the position of the components.
- Remove cement excess before it is completely cured using a suitable instrument (e.g. Le Cron). The components are held in place with light pressure in the process.

Note:

The flow of AGC Cem cement into the space between primary and secondary parts during cementation must be prevented.

- Remove excess material after complete curing using a suitable rotary instrument.

Warning

- Prevent any contact of uncured AGC Cem with the skin / mucous membrane and eyes.
- Uncured AGC Cem may have a slight irritating effect and cause sensitization to methacrylates.
- Customary medical gloves do not provide protection against the sensitizing effect of methacrylates.

Shelf life and storage

- AGC Cem must not be used after the expiration date
- Date of expiration: see note on the syringe
- Storage temperature of AGC Cem: 12–25 °C
- AGC Cem must only be used at room temperature (cooled material is more viscous and takes longer to cure).
- Leave the mixing tip on the syringe after each use. The material will polymerize within the mixing tip, and the mixing tip will therefore serve as a seal.

The material has been developed solely for use in dentistry. Processing should be carried out strictly according to the Instructions for Use. Liability cannot be accepted for damages resulting from failure to observe the Instructions or the stipulated area of application. The user is responsible for testing the products for their suitability and use for any purpose not explicitly stated in the Instructions. Descriptions and data constitute no warranty of attributes and are not binding.

The Material Safety Data Sheet is available online at www.ivoclarvivadent.com

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